

IceBand®

by MD R. Ihrman Patent Pend.

Postoperative cooling bandage

IceBand is a Swedish clinical cold compression bandage specifically developed to relieve postoperative pain and swelling.

IceBand uses cryotherapy/cold therapy to provide an efficient, user friendly and cost-effective method for treating patients after surgery and/or during rehab. Thanks to its unique design, IceBand provides a long-lasting cooling effect for optimal therapeutic benefits, lowering the skin temperature without any risk of frostbites.

IceBand has so far been used to successfully treat more than 1 000 000 patients in Europe. Tests have shown that 93% of the operated patients recommend the usage of IceBand.



Hip/Lower back

- Facilitates rehab
- Reduces pain and swelling
- Cost-effective treatment
- Single-patient use for hygienic safety



IceBand Hip/Lower back



IB 4009

IceBand Hip/Lower back is the newest IceBand model, developed for treating pain and swelling in the hip and lower back.

IceBand Hip/Lower back provides efficient cooling treatment and is suitable for treating postoperative pain and swelling e.g. after hip replacement. The bandage provides cooling and compression for the lower part of the back as well as the hip.

Safe and cost-effective treatment


IceBand has been specially developed as a cryo-therapeutic aid for use after orthopedic surgery. In comparison to other medical devices, IceBand offers a cost-effective approach to treating pain and swelling, enabling orthopedic surgeons to provide efficient post-operative treatment to a large number of patients at a low cost.

Furthermore, IceBand facilitates the rehabilitation and recovery processes, potentially reducing the time spent in the hospital after surgery and optimizing overall hospital care. The compression bandage is designed for single-patient use, significantly minimizing the risk of cross-contamination.

Improved rehabilitation

IceBand should be promptly applied after surgery or injury. When used correctly, it has the potential to enhance postoperative rehabilitation, shorten the recovery period, and diminish the requirement for anesthetics and painkillers. IceBand is suitable for use during physical therapy and rehabilitation training. Additionally, it serves as an excellent aid in the treatment of sporting injuries, such as sprained ankles.

IceBand provides a documented cooling effect for 45-60 minutes. Thanks to the protective non-woven fabric of IceBand, skin temperature will never drop below 50°F (10°C), which rules out any risks for frostbites.



IceBand® is an invention developed in Sweden.
IceBand® is a registered trademark.

IceBand®
by MD R. Ihrman Patent Pend.

Fully MDR-compliant supplier of orthopedic postoperative equipment.

USA: IB Medical Inc., 1643 N. Alpine Road, Suite 104, Box 110, Rockford, IL 61107,
jeff.wagner@icebandmedical.com

Head office: IB Medical AB, Rundelgatan 2, Tranås, Sweden, info@ibmedical.se
Phone: +46 708 14 14 12

www.icebandmedical.com